
Endo-Modell® Arthrodesenagel SK
Modular System

Endo-Model® Knee Fusion Nail SK
Modular System

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CE 0123

WALDEMAR LINK GmbH & Co. KG
Barkhausenweg 10 · 22339 Hamburg, Germany
P.O. Box 63 05 52 · 22315 Hamburg, Germany
Tel.: +49 40 53995-0 · Fax: +49 40 5386929
E-mail: info@linkhh.de · Internet: www.linkhh.de


Endo-Modell® Arthrodesenagel SK

Modular System

Endo-Model® Knee Fusion Nail SK

Modular System

02	Systembeschreibung	02	System Description
04	Indikationen/Kontraindikationen Problemstellung	04	Indications/Contraindications Problems involved
	Implantate		Implants
05	Zementierbare Version	05	Cementable Version
07	Zementfreie Version	07	Cementless Version
08	Sonderanfertigungen	08	Custom-made Implants
	Instrumente		Instruments
12	Instrumente	12	Instruments
	OP-Technik		Surgical Technique
16	Zementierbare Version	16	Cementable Version
24	Zementfreie Version	24	Cementless Version
25	Sonderformen der Distanzüberbrückung	25	Special cases of bridging joint-space defects
27	Fallbeispiel	27	Case History
	Zubehör		Accessories
28	Röntgenschablonen	28	X-ray Templates
	Literatur		Literature
28	Sonderdrucke	28	Special Prints
29	Biomechanische Testung	29	Biomechanical Testing
32	Index	32	Index
	Wichtige Hinweise		Important Information

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The **LINK® Endo-Model® Knee Fusion Nail** consists of a **femoral (1)** and a **tibial (2)** component. A special coupling ensures stable linkage* of the two components. The coupling is achieved with only a small difference in longitudinal travel (about 3 mm) of the Fusion Nail components in relation to one another. After the Fusion Nail components have been inserted in the tibia and femur, the small longitudinal travel facilitates the intraoperative assembly. The two coupling components are semi-circular on one side and have an oblique plane **(3)** on the other. This means that when the components are fitted together by engaging each of them (about 3 mm) in the ring-shaped pocket **(4)** located on the other component, a conical effect is achieved which creates a stable, force-locking connection.

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The connection of both components is finally secured by two **Assembly Screws (5)**. There is one hole in the threaded section of each screw.

A **UHMWPE Lock Bolt (6)** must be placed in each of these holes prior to the final insertion of the screws to secure the screw in position (but not for trial insertion).

The Knee Fusion Nail is produced as **cementable** as well as **cementless version**.

The **cementable version** is made of CoCrMo-alloy or Tilastan® (custom-made for allergy patients). The conical stems are triangular in diameter to provide rotational stability in the cement bed. Femoral stems feature an anatomy adjusted bow. The cementable version of the Knee Fusion Nail includes two UHMWPE centralizers to be screwed upon the end of both stems for correct positioning (see page 06).

The **cementless version** is available as custom-made implant from Tilastan® and its stems are circular in diameter, conical in overall shape and have a fluted surface (see page 07).

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■ Indications/Contraindications - Problems involved

Secondary arthrodesis eg after failed knee arthroplasty or osteosynthesis with loss of function of the extensor apparatus is the only surgical possibility to stabilize the joint and hence maintain function of the lower extremity.

The common procedures to date to perform arthrodesis always depend on bone healing or bony fusion. The variety of methods indicate that this concept very often showed unsatisfactory results. The load-bearing bony substance available to achieve a “natural” arthrodesis through bone healing is considerably reduced and combined with substantial leg shortening, in particular, because of massive loss of bony substance in revision surgery after implant loosening, and extensive debridement in case of infection.

This problem can be resolved by internal stabilization using a Knee Fusion Nail, which can be implanted in the same way as a coupled knee prosthesis. The implant provides primary stability which allows rapid mobilization of the patient. Loss of leg length can also be limited with this concept. In case of infection the cemented Fusion Nail allows specific topical antibiotic application, presenting one of the basic principles of infection therapy.

After repeated exchange surgeries the medullary canal often is enlarged and sclerotic with extensive bone loss in the proximal tibia and the condylar region of the distal femur. Thus, cementing the implant in many cases provides the only possibility for a stable fixation.

Cementless implantation must be viewed under the aspect of secure anchorage - with or without bone grafting.

Contraindications

- Poor general state of health
- Allergies due to (implant) materials.
- Insufficient/inadequate bone mass- or quality which prevents a stable anchor of the prosthesis.
- Distinctive muscular, nerve, vascular or other diseases which put the affected limb at risk.
- Arthrodesis in opposite hip- or knee joint

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Each **cementable Knee Fusion Nail (set)** consists of a femoral and a tibial component, two assembly screws incl. UHMWPE lock bolts and one set (3 ea.) of centralizers. The components will be delivered sterile and in separate packing.

WARNING!

The components of the Knee Fusion Nail SK are not compatible with components of the old version.

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Ø 14 mm montiert, Material: CoCrMo

Femoral components, incl. 1 set of centralizers,

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Tibial components, incl. 1 set of centralizers, set of 2

assembly screws, 4 lock bolts, Ø 14 mm assembled
Material: CoCrMo

Art.-Nr. (Satz) Item no. (Set)	Ausführung Version	Größe Size	Ø A2 mm	B2 mm
	rechts/right links/left			
15-8429/21	X	1	9	170
15-8429/22	X	2	9	210
15-8429/23	X	3	9	250

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Please note:

Cementable Knee Fusion Nails of *Tifootan®* are no standard implants and are available only on request for patients with nickel allergy.

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

The cementless LINK® Endo-Modell® Knee Fusion Nail left and right version from Tilastan® is only available as custom-made Implant.


The production of a cementless Fusion Nail requires full scale X-rays in A-P and M-L planes to decide on its shape and size for optimal implant anchorage. Information about the focal distance and use of a X-ray mounting and measuring bar or correctly on femur and tibia placed metal balls for calculating the magnification on the X-rays is imperative.

Order according to Data Sheet: 160/anagel/4 see page 10.

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Each **cementless Knee Fusion Nail (Set)** consists of a femoral and tibial component, two assembly screws incl. UHMWPE lock bolts.

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■ Custom-made Implants


LINK® Knee Fusion Nail

Available on medical prescription for patients known by name.

Data Sheet: 160/anagel/3

To order please use data sheet. Please copy or request.

 WALDEMAR LINK GmbH & Co. KG Berkhausenweg 10, D-22339 Hamburg Postfach 63 05 52, D-22315 Hamburg Telefon (040) 5 39 95 -0, Telex 2 174 591 link d Telefax (040) 5 38 69 29	LINK® Arthrodesenagel, zementierbar - CoCrMo - LINK® Knee Fusion Nail, cementable - CoCrMo -
---	--

Zentriersterne austauschbar /
centralizers interchangeable

Datum / Date _____

Krankenhaus / Hospital _____

Patient / Patient _____

Termin / Date of Operation _____


Sonstiges/Remarks :

Arthrodesenagel		
Knee Fusion Nail		
Art.-Nr.	Seite	Bestellung
Item No.	Side	Order
	rechts/right	
	links/left	
bestehend aus/consisting of : 1x Femurkomponente/femoral component 1x Tibiakomponente/tibial component 2x Befestigungsschraube/Fixation Screw 2x Kunststoffstopper/HDPE Lock Bolt		

Ein- und Ausschlaginstrument, leihweise		
inserter instrument, on hire		
Art.-Nr.	Benennung	Menge
Item No.	Name	Quantity
15-8450/03	Einführungsinstr. f. fem. Komp.	
	inserter f. femoral component	
15-8450/14	Einführungsinstr. f. tib. Komp.	
	inserter f. tibial component	
15-8450/05	Montageschraube	
	screw for inserter	
15-2534/15	Gewindestab mit Griff	
	threaded rod with handle	
15-8450/06	Schraubendreher	
	screw driver	
131-385	Schraubendreher	
	screw driver	

Zentriersterne		
centralizers		
Art.-Nr.	Menge	Bestellung
Item No.	quantity	Order
15-2975/01	1 Satz / 1 set	
15-2975/01	2 Satz / 2 set	

Material
Materials

CoCrMo

- unmaßstäblich
no scale

- alle Maße in mm
all measurements in mm

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Insertor-Extractor

for femoral component of Knee Fusion Nail with thread, 90 mm

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Insertor-Extractor

for tibial component of Knee Fusion Nail without thread, 90 mm

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Temporary screw for insertor

(only for coupling insertor to nail component – not for connecting femoral and tibial components!)

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Set (3 ea.) **metal trial centralizers**, Ø 12, 14, 16 mm, screwed onto cementable standard Knee Fusion Nails for determine the size of UHMWPE centralizers to be implanted.

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Tibial component of Knee Fusion Nail with inserter and two screws.

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
The threaded rod with handle is screwed into the hole of the inserter and the Knee Fusion Nail is then driven into the medullary canal.

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The following surgical technique refers to an idealized case with minimal amount of bone loss in order to describe all surgical steps conceivable in such a situation. Part of these steps are not applied, of course, if a wide joint gap has to be bridged over in a revision situation, due to heavy loss of bony substance.

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

■ Cementable Version

After opening the knee joint, and the leg in extension, opposite areas of the femur and tibia in correct rotational position are marked. This is best achieved using a colored marker pen or by making a notch on the bone using an osteotome or a rongeur.

After marking the entrance point with a bone awl the tibial medullary canal is opened with an Ø 8 mm twist drill.

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

The entrance point of the femoral medullary canal is marked with a bone awl as well and opened with the Ø 8 mm twist drill.

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Both femoral and tibial joint (shaft) surfaces are resected plane-parallel in such a way as to achieve sufficient surface contact of vital bone tissue between femur and tibia.

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Using long **flexible intramedullary-** or **ball reamers**, tibia and femur are prepared in 1 mm increments to receive the nail components. The minimum diameter of the final reamer corresponds to the diameter of the tip of the nail **A1** respectively **A2** (see opposite drawing and table on page 5). If the medullary canal is wide UHMWPE centralizers of dia. 12, 14 or 16 mm are available to position the nail component in the middle of the canal. Performing the trial implantation trial metal centralizers should be applied.

The depth of the bore is calculated adding to dimension **B1** respectively **B2** (see opposite drawing) plus 15 mm when using a centralizer.

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
The entrance point to take the central section of the Fusion Nail is enlarged further circular around the prebored holes in the tibia and femur, using a gouge and/or a reamer. The diameter of these openings corresponds to the size of the flange on the central section of the Fusion Nail (see drawing on page 6). The required depth of these openings in the femur and tibia is 25 mm each. The central section of the nail (dimension D on the nail dimension sheet), with a total length of 50 mm, is thus to be inserted with one half in the femur and the other in the tibia.

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Drawing referring to dimensions
A, B, C, D on page 6

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Ventrally a 30 x 30 mm fenestration in axial direction has to be prepared with an oscillating saw.

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Using the **insertor** attached to the **threaded rod and handle** the tibial and femoral components of the Fusion Nail are introduced into the prepared cavities. Femur and tibia are then positioned on each other while, at the same time, fitting both components together on a trial basis. If it should prove difficult to insert either of the components it may be necessary to enlarge the medullary canal using an **intramedullary reamer**, a **ball reamer** or a **curette**.

After the trial repositioning and before removing the Fusion Nail components marks are placed on the implants as well as on the bone to determine the position of the implants for the final insertion. The insertion depth is also marked on the implants.

■ Surgical Technique

This ventral opening is necessary in order to be able to place the **assembly screws** ventrally after the nail components have been driven home completely. The resected osseous lamellae are then used to later close the ventral fenestration.

If femur and tibia cannot be joined together the remaining gap is bridged only by the Fusion Nail itself or filled in addition with bone graft. A ventral fenestration does not apply if the gap measures 50 mm or more or its height has to be reduced according to the smaller distance (see special cases on page 25).


The nail components are inserted into the prepared bone cement in accordance with the markings which have been made on them (whether femoral or tibial component first, depends on the individual situation).


Note that the inserter must be removed from the nail component 2 to 3 cm before the latter has reached its final position. The Fusion Nail is then pressed by hand into its final position.

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After cement hardening both components of the Fusion Nail are joined together and secured with two **assembly screws** with the UHMWPE lock bolts placed in the threaded shaft before. If bone contact between femur and tibia can be achieved the bone surfaces must be free of cement. To facilitate bony union the contact surfaces should be freshened. If appropriate, additional bone material can be introduced between the contact surfaces.

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The gap to be bridged has to be filled preferably with cement if bony contact between femur and tibia cannot be achieved.

Before joining the components anchoring holes for the cement used to fill the gap (abt. 5-6 mm wide and abt. 10 mm deep) should be prepared in the contact surface areas to improve the bond between femur and tibia.

It is necessary for the surgeon to get familiar with the Fusion Nail before the surgery. The Knee Fusion Nail is delivered sterile.

NOTE:

The UHMWPE lock bolts for the assembly screws are delivered sterile within a transparent pack: They are visible once the external packaging is removed. The transparent packaging of the UHMWPE lock bolts must be opened, however, only during the operation under sterile conditions.

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
■ Cementless Version

The surgical technique for the **cementless** and **cementable** versions is identical except for the impaction of the Fusion Nail components.

Both sections of the Fusion Nail components are driven in position carefully until the desired position between femur and tibia is obtained by means of repeated trial reductions.

Once the **Fusion Nail** components have been driven into the desired position the securing assembly screws are screwed in. Before finally turning in the assembly screws the UHMWPE lock bolts must be inserted into the assembly screws to secure their fixation. The assembly screws must then be tightly screwed in. The bone block earlier removed is used to fill the ventral window.

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
■ Special cases of bridging joint-space defects


A remaining space between femur and tibia, bridged by the Fusion Nail and possibly filled with bone graft reduces the necessary depth of the implant bed in femur and tibia, in case of a 50 mm distance, altogether.

Example:

The central part D of the Fusion Nail is 50 mm. If a gap of 20 mm between femur and tibia is to be bridged the depth of the implant bed in femur and tibia comes up to 15 mm each, altogether 30 mm. According to the situation it may be advisable, however, to plan the total depth of 30 mm for only the femur or tibia, or something in between.

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■ Case History

54 years old female patient,
suffering 28 years of chronic polyarthritis;
multiple pre-operations,
including arthroplasties in various joints.

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3/93 patient arrived with a centrally located crater of approximately
one inch diameter within the knee joint. Total painful loss of func-
tion in the extremity.

Infection with staphylococcus aureus and pseudomonas aerugi-
nosa.

Soft tissue defects and substantial loss of bone stock following
removal of prosthesis.

3/93 Revision mit Implantation des LINK® Endo-Modell®
Arthrodesenagels.

3/93 revision with implantation of LINK® Endo-Model® Fusion Nail.

- primäre Stabilität
primary stability
- primäre Wundheilung
primary wound healing
- Mobilisation unter Vollbelastung
mobilization under full load
- Entlassung 44 Tage postoperativ schmerzfrei, voll mobilisiert
discharge 44 days postoperatively, painless, fully mobilized
- 1999 weiterhin beschwerdefrei
follow-up 1999. Patient without complains.

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Cemented modular intramedullary nail in failed knee
arthroplasty-a report of 2 cases
Acta Orthop Scand 1999; 70 (3): 305-307

Wodtke JFK
Endo-Model® Knee Fusion Nail
Waldemar Link GmbH & Co. KG, (2000) Endo-Klinik Hamburg, Germany

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

■ Biomechanical Testing

Dynamic study of the clamping assembly of a Knee Fusion Nail

(Dip.-Ing. Georg Duda, Technical University Hamburg/ Harburg, Biomechanics Department, Prof. Dr. E. Schneider, 02/94)

Extract from the Study

Questioning

How does the coupling device of the Knee Fusion Nail behave under alternating loads, with particular fatigue and surface wear. What effect do screw starting torque, surface condition, and manufacturing and fit tolerances have on the dynamic behavior of the Knee Fusion Nail.

Material and Methods

The study material is divided into two comparable groups of nails that are equal in size. The „high cycle“ fatigue test is expected to produce evidence with regard to the fatigue strength and form strength of the coupling elements. The test involves subjecting the nail to estimated average dynamic loading of a least 1 million-load cycles. The aim of the „high cycle“ fatigue test is to bring about artificial aging of the nail and from this to draw conclusions with regard to an optimal design.

In the „low cycle“ fatigue test, on the other hand, only a limited number of cycles is applied at maximum loading in order to learn about the failure mechanisms and weak points under extreme loads.

Study Material

WALDEMAR LINK GmbH & Co. KG provided eight Knee Fusion Nails for use in the study. The samples were in accordance with the dimensions of standard commercially available devices including the geometry and diameter of the stems.

Test Machine and Data Processing

The study was carried out on a servo-hydraulic Bionix test machine manufactured by MTS of Minneapolis, USA. In all of the tests the assembly screws were tightened at a torque of 11 N-m. The data were evaluated in Microsoft Excel.

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
Abb./Fig. 3.2a


Abb./Fig. 3.2b

It is evident that under the loads applied in this test the wear marks are minimal and that they have no impact on the functional performance of the Knee Fusion Nail.

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Evaluation of the test results in the compression test shows that failure is unlikely in this loading mode.

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The results of the torsion test show that the dynamic behavior of the Knee Fusion Nail does not change after a relatively short run-in period.

„Low cycle“ fatigue test

4-point deflection, Figs. 3.11, 3.12

No failure of the nail connection occurred during 4-point deflection.


Abb./Fig. 3.12

Torsion test

The Modular Knee Fusion Nail exhibits exceptional torsional rigidity due to the geometry of the coupling elements. Even under extremely high loading there was no failure of the coupling elements. Implant failure should be unlikely even under combined flexional and torsional loading.

Conclusion

The Knee Fusion Nail met all test standards without problems. Only during bending parallel to the division plane were there minimal movements between the segments. The Knee Fusion Nail exhibits great axial and torsional rigidity.

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Please note the following regarding the use of our implants:

- 1. Choosing the right implant is extremely important.**
The size and shape of the human bone determine the size and shape of the implant and also limit the load capacity. Implants are not designed to withstand unlimited physical stress. Demands should not exceed normal functional loads.
- 2. Correct handling of the implant is exceedingly important.**
Under no circumstances should the shape of a finished implant be altered, as this shortens its life span.
Our implants must not be combined with implants from other manufacturers.

The instruments indicated in the Surgical Technique must be used to ensure safe implantation of the components.
- 3. Implants must not be reused.**
Implants are supplied sterile and are intended for single use only. Used implants must not be reused.
- 4. After-treatment is also very important.**
The patient must be informed of the limitations of the implant. The load capacity of an implant cannot compare with that of a healthy bone!
- 5. Unless otherwise indicated, implants are supplied in sterile packaging.**
Note the following conditions for storage of packaged implants:
 - Avoid extreme or sudden changes in temperature.
 - Sterile implants in their original, intact protective packaging may be stored in permanent buildings up until the "Use by" date indicated on the packaging.
 - They must not be exposed to frost, dampness and direct sunlight or mechanical damage.
 - Implants may be stored in their original packaging for up to 5 years after the date of manufacture. The "Use by" date is indicated on the product label.
 - Do not use an implant if the packaging is damaged.
- 6. Traceability is important.**
Please use the documentation stickers provided to ensure traceability.
- 7. Further information** on the material composition is available on request from the manufacturer.

Follow the instructions for use!

WALDEMAR LINK GmbH & Co. KG, Hamburg.

All content in this catalogue, including text, pictures and data, is copyright-protected. Every instance of use not permitted by the German Copyright Act is subject to our prior consent. In particular, this applies to the reproduction, editing, translation, saving, processing or passing on of content stored in databases or other electronic media and systems. The information in this catalogue is solely intended to describe the products and does not constitute a guarantee.

The Surgical Technique described has been written to the best of our knowledge and belief, but it does not relieve the surgeon of his responsibility to duly consider the particularities of each individual case.

Unless otherwise indicated, all instruments are made of surgical stainless steel.

WALDEMAR LINK GmbH & Co. KG
Barkhausenweg 10 · 22339 Hamburg, Germany
P.O. Box 63 05 52 · 22315 Hamburg, Germany
Tel.: +49 40 53995-0 · Fax: +49 40 5386929
E-mail: info@linkhh.de · Internet: www.linkhh.de