

Endo-Model[®] Protésis de Rodilla Rotacional y Charnela

Endo-Model[®]
Rotational and Hinge Knee
Prostheses


Técnica Quirúrgica
Surgical Technique

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To initially open the femoral canal, the bone awl is placed at the lowest point (red) of the trochlear groove in the patellofemoral joint. The point lies at the tip of the intercondylar fossa.

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
The femoral canal is opened with an Ø 8-mm twist drill at the point previously marked.

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
The femoral canal is carefully widened with a ball reamer.

15-1133/02E	Fresas esféricas 250 mm - Ø 10 mm
15-1133/03E	Fresas esféricas 250 mm - Ø 12 mm
15-1133/04E	Fresas esféricas 250 mm - Ø 14 mm
15-1133/05E	Fresas esféricas 250 mm - Ø 16 mm
15-1133/06E	Fresas esféricas 250 mm - Ø 18 mm

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

The femur rasp is used to adjust the opening made with the twist drill and ball reamer to fit the shape of the femoral stem.

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The femoral saw guide corresponds in shape to the portion of the femoral component that will lie within the bone (A).


To keep it centered within the canal, the tip of the stem of the femoral saw guide is equipped with a metal centralizer available in different sizes corresponding to the diameter of the femoral canal. **The centralizer serves as a placeholder, simulating the polyethylene centralizer on the implant (B).**

The femoral saw guide is advanced until its box is in contact with the condyles. The rotational adjustment is correct when the anterior groove of the resection box matches the trochlear groove when viewed

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

Anterior and bilateral femoral osteotomies are performed with the oscillating saw along the surface of the box (D).

The resected bone at the junction of the condyles and the femoral diaphysis should be carefully osteotomized step by step. This should allow the box of the femoral saw guide to be inserted up to the level of the condyles after the resected bone has been removed with a lambotte osteotome (E).

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
The anterior femoral canal must be adapted to accommodate the slight angle between the anterior surface of the intracondylar box and the stem of the femoral saw guide and the corresponding angle in the prosthesis. To do this, the cancellous bone posterior to the patellofemoral joint is hollowed out in retrograde fashion with the oscillating saw.

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
Hollowing out the bone as described creates a step-off in the cancellous bone proximal to the articular surface. This step-off must be removed with the reamer or rasp. If the step-off is not removed, the anterior surface of the box of the femoral saw guide will not fit tightly against the anterior bone. This will create a gap anteriorly between the physiologic trochlear groove and its prosthetic continuation.

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The femoral canal is intentionally first opened far posteriorly (Fig. F, red point). The inserted femoral guide must then be shifted anteriorly (Fig. F, blue point) until the stem of the femoral saw guide (locked in extension) is aligned with the distal femur (Fig. G, blue position of femoral saw guide), so as to prevent hyperextension.

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
If proper axial alignment of the femoral saw guide cannot be achieved after steps F3 – F4 have been performed, the femoral saw guide will have to be shifted further anteriorly. This requires resection of more bone from the intracondylar portion of the trochlear groove. A lambotte osteotome impacted into the bone can be used to define the extent of the osteotomy. The impacted osteotome also guides the saw blade until the saw has cut a groove in the bone. It is important to perform the resection gradually, verifying proper alignment of the axis of the femoral saw guide with the axis of the femur.

After this additional resection, the steps of hollowing out the cancellous bone proximal to the trochlear groove and removing cancellous bone from the anterior portion of the femoral canal proximal to the joint must be repeated.

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
The box of the femoral saw guide is advanced into the bone until the contour of the anterior groove of the saw guide is congruent with the trochlear groove in the patellofemoral joint. The femoral saw guide is fixed with 2 fixation pins. The condyles are then adapted to match the curved surface of the resection box.

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
The patellar glide resection guide is attached to the saw guide and the anterior surface of the distal femur is resected to fit the shape of the implant.

15-2622	Raspa, 350 mm Rasp, 350 mm
---------	----------------------------

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
The remaining anterior edge is rounded off with a rasp. Residual cartilage should be removed.

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In sclerotic bone, it is a good idea to achieve a slight curvature of the inner edges of the condyles to match the interior contour of the implant between box and flanges. Soft cancellous bone will conform to the interior contour of the flanges; sclerotic bone must be shaped to fit. Where the medial femoral condyle is sclerotic, cement retention holes are drilled to improve implant fixation.

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
A cylindrical bone plug is inserted into the proximal femur to seal the femoral canal. The bone plug is advanced to a depth corresponding to the length of the stem of the femoral component. This bone plug helps to achieve hemostasis in the proximal medullary canal. It seals the canal when the cement is later injected during implantation of the femoral component.

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
The trial prosthesis is inserted into the femoral canal to check the resection. Once the femoral preparation is complete, the femoral component can already be cemented in at this stage. This avoids accidental fracture of the femoral condyles during the operation.

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B Rotational and Hinge Knee Prostheses without Patella Flange (Fig. F7-b to F9-b)

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

The box of the femoral saw guide is advanced into the bone until the contour of the anterior groove of the femoral saw guide is congruent with the trochlear groove in the patellofemoral joint. The condyles are then adapted to match the curved surface of the resection box.

In sclerotic bone, it is a good idea to achieve a slight sagittal curvature of the inner edges of the condyles to match the interior contour of the implant between box and flanges. Soft cancellous bone will conform to the interior contour of the flanges; sclerotic bone must be shaped to fit. In cases where the medial femoral condyle is sclerotic, cement retention holes are drilled to improve fixation of the implant.

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A cylindrical bone plug is inserted into the proximal femur to seal the femoral canal. The bone plug is advanced to a depth corresponding to the length of the stem of the femoral component. This bone helps to achieve hemostasis in the proximal medullary canal. It seals the canal when the cement is later injected during implantation of the femoral component.

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The trial femoral prosthesis is inserted into the femoral canal to check the resection. The articular surface of the trial prosthesis must lie flush with the articular surface of the trochlear groove, without any gap between implant and bone.

Once the femoral preparation is complete, the femoral component can already be cemented in at this stage. This avoids accidental fracture of the femoral condyles later on during the operation.

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

Next the tibia is exposed. The hand on the opposite side to the approach is placed in the popliteal fossa with the thumb abducted, while the other hand grasps the ankle anteriorly. Applying traction and shear force with the hands with slight external rotation exposes the proximal tibia.

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
The tibial canal is marked with a broach at the junction between the anterior and middle thirds of the sagittal diameter of the tibial surface, immediately anterior to the intercondylar eminence. The point described lies over the center of the medullary canal of the tibia. An 8-mm twist drill (A) is used to open the tibial canal. The tibial canal is then carefully enlarged distally with a ball reamer (B). To seal off the distal tibial canal, autologous cancellous bone is pressed deep into the tibial canal with an impactor, down to a point inferior to the tip of the tibial component stem (C).

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The intercondylar eminencia is removed down to the level of the intact tibial plateau to create a reference plane for measuring the bone to be resected.

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
The tibial saw guide is inserted into the tibial canal, and aligned with the axis of the tibia in the median lateral plane and coronal plane. The height of the tibial saw guide is adjusted such that 10 mm can be resected. This resection level corresponds to the thickness of the thinnest portion of the implant's polyethylene plateau (D).

After the tibial saw guide is removed, the osteotomy is deepened. The resected bone is mobilized and sharply dissected off the posterior capsule with an osteotome (E).

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
The alignment guide is removed from the saw guide to determine rotational alignment.

The mark on the metal plateau of the tibial saw guide and the horizontal part of its anterior alignment rod will indicate an area between the middle of the tibial tuberosity and its medial margin (F). The second digital ray of the foot should exhibit external rotation between 0° and 10° (G).

The rotational alignment of the tibial component corresponding to the mark on the metal plateau is best marked on the anterior cortex of the proximal tibia by a small notch made with a narrow bone rongeur (F).

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
The tibial broach, which is shaped like the tibial stem of the implant, is used to shape the tibial canal. The anterior pin on the rasp is aligned with the notch made in the proximal tibia. The tibia reamer is successively driven in and withdrawn until it can be inserted down to the level of the resection plane.

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
Optionally trial femoral and tibial prostheses can be inserted and coupled for trial reduction.

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
It is now important to verify correct sagittal and rotational alignment of the implant before cementing the implant components in place:

- A** hyperextension
- B** impaired extension
- C** external hyperrotation
- D** internal hyperrotation


must be corrected.

Figures **E** and **F** show correct axial alignment of the seated implant.

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
Once it has been verified that the implant is correctly positioned, the components are individually cemented in place. It is recommended that the femoral component normally be implanted with 80 g of bone cement. Excess cement at the sides is removed.


15-2538	Impactor para Componentes Femorales - pequeño / mediano
15-2537/02	Impactor para Componentes Femorales - grande

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The plastic strip that protects the joint mechanism from the entry of cement is removed. The strip should be pulled posteriorly.

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
Before the tibial component is cemented with at least 40 g of cement, the tibial canal is sealed with a bone plug. Note the marked rotational position when aligning the implant.

Caution: The tibial component may only be cemented without polyethylenen plateau removed where the trial screw has been inserted to a maximum depth.

To remove the polyethylene plateau, the trial screw is removed and the plateau taken off the metal tibial tray using the tibial plateau introducer. Subsequently the trial screw is reintroduced until stop in the metal plateau. **This is necessary to prevent cement from entering the threaded hole in the metal plateau.**

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Before the cemented components are fitted together, the polyethylene plateau still has to be removed from the tibial metal tray. With the knee flexed, the femoral component is inserted onto the pin of the tibial component.

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
To insert the polyethylene tibial plateau, the femoral component is lifted slightly. The tibial plateau is then inserted from an interior direction between the proximal and distal component of the prosthesis (G).

Care should be taken to ensure that the chamber of the plastic plateau merges with the flange of the femoral component, and that the dovetail recess on the underside of the plastic plateau snaps into the marginal groove on the metal tibial tray (H).

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The polyethylene plateau is secured in place on the metal tibial tray with the self-locking fixation screw.

Caution:


The self-locking fixation screw must only be used during the final assembly of the plateau. Loosening the fixation screw destroys the screw retention system in the polyethylene plateau, and a new plateau must then be inserted.

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The implanted Rotational Knee Prosthesis should allow up to 90° of flexion, depending on the soft tissue structures in extension, a resilient extension impairment of approximately 5° is optimal. This extension impairment helps to ensure secure contact between the two prosthesis components.

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
**X-ray Templates for
Endo-Model® Total Knee Prostheses
(Rotational and Hinged version)**

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
Information on Implants and Instruments
Catalog available on request

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Important Information about our Orthopaedic Implants

More detailed information about our implants can be supplied upon request.

Materials used for our orthopaedic implants:

- CoCrMo alloy, ISO 5832-4/ASTM F75 or ISO 5832-12
- TiAl6V4 alloy, ISO 5832-3/ASTM F136 or ASTM F1108
- Unalloyed titanium, ISO 5832-2/ASTM F67
- Stainless steel, ISO 5832-1/ASTM F138 / ASTM F139
- CoCrNiMoFe, ISO 5832-7/ASTM F1058
- Al₂O₃ (aluminium oxide ceramic), ISO 6474
- UHMWPE, ISO 5834-2/ASTM F648
- Calcium phosphate coating, ASTM F1609
- Hydroxyapatite coating, ASTM F1185

Please note the following regarding the use of our implants:

1. Choosing the correct implant is extremely important. The size and shape of human bones restrict the size and shape of the implant. This also restricts the implant's ability to withstand stress. Implants are not designed to withstand unlimited body loads. The physical demand on implants should not surpass normal functional stresses.

2. It is also very important to treat implants correctly. Altering the shape of a finished implant reduces its lifespan. The following must be actively avoided: bending the implant sharply, bending it back, or making grooves or scratches in it. Our implants may not be combined with implants from other manufacturers.

3. Implants must not be re-used. Even if a used implant looks undamaged externally, it must be assumed that the material has become fatigued internally.

4. Patient rehabilitation is also very important. The patient must learn about the limitations of the implant. An implant cannot compare with healthy bone in the ability to withstand stress!

5. Unless otherwise stated, implants are shipped in sterile packaging. The following are important for storage of packaged implants:

- Avoid extreme or sudden changes in temperature.
- We recommend a storage temperature of 18 –20 °C with 50 – 70% humidity.
- Avoid direct sunlight.
- Protect implants from damp and from mechanical damage.
- Implants may be stored for up to 5 years in their original packaging. The 'Use by' date is on the product label.
- Do not use an implant if its packaging has been damaged.

6. The documentation trail is important. Please use the documentation sticker included in the packaging for this purpose.

WALDEMAR LINK GmbH & Co KG, Hamburg

All materials in this catalogue, including text, pictures and data, are protected by the relevant copyright provisions. Every use of this material that is not allowed by the German Copyright Act is subject to our prior consent. This especially applies for reproduction, alteration, translation and or publication of electronically or otherwise saved, modified and/or published contents. All information contained in this catalogue is solely intended for the description of products and shall under no means establish any kind of warranty.

All instruments, unless otherwise noted, are made of Stainless Steel.

WALDEMAR LINK GmbH & Co. KG

Barkhausenweg 10 · D-22339 Hamburgo

Apartado Postal 63 05 52 · D-22315 Hamburgo

Teléfono : +49 (0)40 5 39 95-0 · Fax: +49 (0)40 5 38 69 29

e-mail: info@linkhh.de · Internet: www.linkhh.de